

```
>> OSA202C_TestCode
```

```
osa =
```

```
LibDeviceInterface with properties:
```

```
    AutoStoreOptions: [1x1 ThorlabsOSAWrapper.AutomaticFileStoreOptions]  
InstrumentDescription: [1x1 ThorlabsOSAWrapper.DeviceProperties]  
    AcquisitionSettings: [1x1 ThorlabsOSAWrapper.DeviceAcquisitionSettings]  
        MotorRunning: 1  
    IndexOfSpectrometer: 0
```

```
Connected to OSA
```

```
Setting Sensitivity Mode
```

```
ans =
```

```
OperationResult with properties:
```

```
    Failure: 0  
    Success: 1  
ErrorMessage: []  
    ReturnCode: 0
```

```
Setting Resolution Mode
```

```
ans =
```

```
OperationResult with properties:
```

```
    Failure: 0  
    Success: 1  
ErrorMessage: []  
    ReturnCode: 0
```

```
Setting Setting Auto Gain
```

```
Is auto gain on?: true
```

```
Acquiring Spectrum.
```

```
ans =
```

```
OperationResult with properties:
```

```
    Failure: 0  
    Success: 1  
ErrorMessage: []  
    ReturnCode: 0
```

```
Spectrum Acquired
```

```
ans =
```

OperationResult with properties:

Failure: 0
Success: 1
ErrorMessage: []
ReturnCode: 0

Connection to OSA closed

Length of measured data: 88077

- If non-zero, data is probably present

Attempt to create SpectrumStruct

Error using OSA202C_TestCode (line 43)

Message: The namespace for the type 'spectrum_t*' does not exist. MATLAB requires all
NET types to
be within a namespace.

Source: dotnetcli

HelpLink:

>> disp('Try Other stuff')

osa.GetLastSpectrum() % call without provided struct

Try Other stuff

No method 'GetLastSpectrum' with matching signature found for class
'ThorlabsOSAWrapper.LibDeviceInterface'.

>> osa.GetLastSpectrum(ThorlabsOSAWrapper.SpectrumStruct(osa.GetLastSpectrumLength())) %
creat struct in-line

Message: The namespace for the type 'spectrum_t*' does not exist. MATLAB requires all
NET types to
be within a namespace.

Source: dotnetcli

HelpLink:

>> asm2 = NET.addAssembly('C:\Program Files\Thorlabs\Thorlabs OSA\FTSLib.dll');

Message: The given assembly name or codebase was invalid. (Exception from HRESULT
0x80131047)

Source: mscorlib

HelpLink:

>>